Viscan AB-101 KULLANICI KILAVUZU

DIKKAT

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Estimated waist circumference


■ Reliability of estimated waist circumference measurement

The AB-101 uses the high correlation between waist circumference at the navel as measured using a tape measure (abdominal circumference) and the breadth of the abdomen to calculate the waist circumference. The regression equation is formulated based on data from approximately 600 Japanese subjects of a variety of physiques from infants to athletes, healthy adults and people with spinal injuries, etc.

The estimated waist circumference measurement achieves high levels of repeatability and excellent reliability by utilising a positioning laser and multiple non-contact sensors.

*The estimated waist circumference is a value for the waist circumference at the navel (abdominal circumference). This unit can be used to track changes in subjects for whom measuring waist circumference has previously been difficult.

Correlation between waist circumference as measured by tape measure and estimated waist circumference according to the AB101.

* This scatter graph is based on 120 healthy males and females.

Internal organ fat levels

■ Standards for internal organ fat levels

Internal organ fat levels are best used to track changes over the long term for effective health management.

Internal Organ Fat	1.0 – 5.0	5.5 – 9.5	10.0 – 12.0	12.5 – 14.5	15.0 – 17.0	17.5 and more
Bar number	1	2	3	4	5	6
Determination	Average		Slightly excessive		Excessive	
Interpretation	No need for concern at present. Continue with a balanced diet and an appropriate amount of exercise.		Ensure that an appropriate amount of exercise is done and calorie intake is limited to reduce weight to an appropriate level.		Weight loss is required through actively engaging in exercise and food restrictions. Please see a medical practitioner for a medical diagnosis.	

* Medical opinion states that when the actual value is more than 100cm_ according to X-ray and CT imaging there is a high probability of resultant life-style related disease. A level of 10.0 with this unit generally correlates to an internal organ fat area of 100cm_.

Internal organ fat levels measured by this unit.

- It is possible for subjects with low central fat percentages to have high internal organ fat levels.
- Please see a medical practitioner for a medical diagnosis.

Central fat percentage

What is the central fat percentage?

True obesity is determined not by weight but by central fat percentage. Not only is body fat a source of energy for activity, it also maintains body temperature, protects the body from external impact, provides moisture for our skin, creates a smooth body line and ensures normal hormone activity, among other functions. An excessively high percentage of body fat, however, causes hypertension, hyperlipidemia, diabetes and other life-style related conditions.

There is said to be a particularly strong link between life-style related conditions and body fat in the trunk (central fat ratio). This unit analyses data based on dual energy x-ray absorptiometry (DXA) which reveals the composition of individual body segments, and determines the central fat percentage using this analysis as a benchmark. See below for central fat ratio determinations.

* The 'central body' referred to here consists of the head, neck, chest, abdomen and lower back.

Central fat determination

[Males]



Central Fat (%)	Standard Deviation	Bar Number	Determination
<7.3	up to -2SD	2	Low
$7.3 \leq - < 13.9$	-2SD - -1SD	3	
$13.9 \leq - < 17.2$	-1SD - -1/2SD	4	
$17.2 \leq - < 23.8$	-1/2SD - +1/2SD	5	Average
$23.8 \leq - < 27.1$	+1/2SD - +1SD	6	
$27.1 \leq - < 33.7$	+1SD - +2SD	7	High
$33.7 \leq - < 40.3$	+2SD - +3SD	8	
$40.3 \leq$	+3SD or more	9	

[Females]

Central Fat (%)	Standard Deviation	Bar Number	Determination
<6.0	up to -3SD	1	Low
$6.0 \leq - < 13.7$	-3SD - -2SD	2	
$13.7 \leq - < 21.5$	-2SD - -1SD	3	
$21.5 \leq - < 25.3$	-1SD - -1/2SD	4	Average
$25.3 \leq - < 33.0$	-1/2SD - +1/2SD	5	
$33.0 \leq - < 36.9$	+1/2SD - +1SD	6	High
$36.9 \leq - < 44.6$	+1SD - +2SD	7	
$44.6 \leq - < 52.3$	+2SD - +3SD	8	
$52.3 \leq$	+3SD or more	9	

* Categories are determined according to average values of adult population samples and standard deviations.

Correlation between DXA central fat percentage and abdominal BIA central fat percentage

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TEKNİK ŞARTNAME

Model		AB-101	
Impedance Meter	Measuring method	4 electrode format	
	Measuring frequency	6.25kHz, 50kHz	
	Measuring current	500 μ A	
	Electrode material	Stainless steel	
	Measuring site	Abdomen (at navel)	
	Measuring range	5 – 100 Ω	
	Power source	DC 6V/4 AAA alkaline batteries (included)	
Main Unit	Measuring method	Near infrared reflectance	
	Measuring range	20 – 48cm (correlates to estimated waist circumference of 50 – 130cm)	
	Battery	Nickel-hydrogen rechargeable batteries; output DC 7.2V, 1mAh	
	Power source	AC adapter Input: 100 – 240V Output: 12V, 2A Plug: centre minus	
Laser (for positioning)	Light emitter	Visible light semiconductor laser	
	Laser wavelength	Wavelength: 650 nm (red light)	
	Laser output	Max. 1mW	
	Beam size	At a distance of 200mm: minimum line length 200mm; light width max. 2mm Point \pm 5mm from the centre (point is approx. \varnothing 1.5).	
	Classification	Class 2M	
Display	Estimated waist circumference	50 -130cm (increment: 1cm)	For subjects 18 years and over
	Central fat %	Increment: 0.1% (5.0 – 75%)	
	Internal fat level	1.0 – 59.0 levels (increment: 0.5 level)	
	Central fat % determination	Low – high (with bar indicator)	
	Internal fat determination	6 levels from average – excessive (with bar indicator)	
Settings	Male / Female		
Total Weight	Main unit: approx 2.6kg; Impedance meter: approx. 0.5kg		
External interface	RS232C output (D sub 9 pin connector, female)		
Usage temp. range	5 - 35°C		
Temperature range for storage	-10 - 50°C		

İRTİBAT:

TARTI DIŞ TİC.ve PAZ.LTD.ŞTİ.

Dikilitaş Mahallesi, Karanfil Sokak , Kartal Apt. No.15
34330 Yıldız - Beşiktaş / İstanbul

Tel : 0212 327 10 05

Faks : 0212 327 09 98

www.tarti.com info@tarti.com